SINERGIA DISK Ti

Ti6Al4V-ELI

Ti	Al	V	Altri - Others - Andere - Autres - Otros
90%	6%	4%	< 0.4%

ISO 22674 – Type 4

CE
0546

Caratteristiche tecniche • Technical data • Technische Daten Caractéristiques techniques • Características técnicas

Densità - Density - Dichte - Densité - Densidad	g/cm ³	4,4
Intervallo di fusione - Melting range - Schmelzintervall - Intervalle de fusion - Intervalo de fusión	°C	1605 - 1660
CET - CTE - WAK - CET - CET 25-500 °C - 25-600 °C	10 ⁻⁶ K ⁻¹	9,6
Limite elastico 0,2% - Yield strength 0,2% - 0,2% Dehngrenze - Limite élastique 0,2% - Limite elástico 0,2%	MPa	>760
Allungamento - Elongation - Bruchdehnung - Allongement - Alargamiento	%	>8
Modulo elastico - Modulus of elasticity - Elastizitätsmodul - Module d'élasticité - Módulo de elasticidad	GPa	115
Durezza Vickers - Vickers hardness - Vickershärte - Dureté Vickers - Dureza Vickers	HV	350

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SINERGIA DISK Ti are disks made in a Titanium based alloy of excellent quality (Ti6Al4V-ELI) for carrying out dental prostheses by CAD/CAM milling centers.

Ti6Al4V-ELI is a Titanium Grade 5 with highest purity ("Extra-Low Interstitial") and is used since decades in the biomedical field thanks to its excellent biocompatibility, corrosion resistance and physical-mechanical properties.

INDICATIONS

Suitable for abutments, single crowns and long-span bridges on natural dies or implants:

- Anterior: not longer than 3-units between the pillars, connection section > 6 mm²
- Posterior: not longer than 3-units between the pillars, connection section > 9 mm²

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Finishing: use carbide disks (or other proper instruments) to separate the milled framework from the disk. For the finishing use only suitable burs in tungsten carbide with cross-cut and work in the same direction. The milled frameworks can be pre-polished with silicone rubbers and polished with normal polishing pastes and soft brushes (goat hair).

NOTE! During the separation from the disk and the finishing phase, do not overheat the metal frame by avoiding too high speed and pressure. Often cool the piece in water. Never use hydrofluoric acid.

Aesthetic finishing: sandblast with aluminum dioxide 125 - 250 µm at 2-3 bar and further clean with steam. Use composite or ceramic suitable with Titanium. The CTE of the ceramic should be compatible with Titanium (25-500°C = $9,6 \times 10^{-6} \text{ K}^{-1}$), the firing temperature should not be higher than 800°C. Follow the instructions for use of the ceramic manufacturer.

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